

Step out with confidence

The all-new Class 3 **ST#5D** scooter has a modern, robust design that gives you the latest in scooter styling and performance. With its bespoke alloy wheels, sleek paintwork and touch screen controls, it's a real trendsetter. Despite being one of the larger scooters in the range, it's compact and easy to manoeuvre with a tight turning circle.

Modern touch screen controls

For the ultimate modern look and feel, the **ST#5D** has an easy to use LED touch screen control panel. Simply press the screen to change your speed, use your horn or activate your indicators.

Advanced suspension for maximum comfort

The **ST#5D** isn't just stylish – it offers outstanding performance too. Using proven technology from the motorcycle industry, this revolutionary scooter has independent left and right front suspension and single centred mono-shock suspension to the rear. This, combined with a high ground clearance gives you an incredibly smooth ride, even on rougher terrain.

The **ST#5D** has the same standard features as the **ST#1**, plus the following:

Additional Features

- Increased range and weight capacity
- Suitable for use on public highways
- Luxurious fully and easily adjustable seat with a headrest
- Two rear view mirrors
- LED indicators, rear lights and brake lights
- Tiller adjustment
- Delta handlebar
- Touch screen control panel
- Thumb lever driving controls
- High / low speed selector
- Advanced suspension
- Emergency brake
- Bespoke alloy wheels with pneumatic tyres
- Front shopping basket
- Increased ground clearance
- Power / self-diagnostic warning light
- Rhino 2 control system
- Tie down points for non-occupant transportation

Strider

ST#5^D

Strength with stability



28ST
180KG



8MPH
12KPH



26.1MILES
42KM



12 MONTH
WARRANTY



ANTHRACITE
GREY AVAILABLE
PAINT
OPTION