

A New Age In Scooter Design...

...The Next Generation Drive Scooter.



www.cobrascooter.co.uk

drive
www.drivemedical.co.uk

The Cobra is a revolutionary scooter in terms of design, aesthetics and styling. It has been aerodynamically designed to give an outstandingly stable drive and maximise battery consumption. It's a remarkable scooter that makes a huge statement and sets a new benchmark with its striking dynamic design.



Revolutionary design, styling & performance

- Top speed of 8mph*
- 182kg (28st) weight limit
- Delta bar
- Excellent turning circle
- Rotating and sliding seat
- 3 colours available
- 30 mile range*
- 14" diameter wheels
- 2 x 70ah batteries
- Twin rear view mirrors
- Road & pavement legal
- Height adjustable headrest
- Width and angle adjustable armrests
- All round suspension
- LED lighting
- High Level Indicators
- Brake system / electromagnetic & emergency brakes
- 'Blade' alloy wheels
- Anti tip wheels
- Low profile tyres

Colour Options



WHITE



GRAPHITE



BLACK


SPECIFICATION



STATE OF THE ART CONTROLS



'BLADE' ALLOY WHEELS




LED LIGHTING

Code: COBRA
Description: COBRA




LENGTH
155cm (61")



HEIGHT
135cm (53")



WIDTH
70.5cm (27.6")



MAX CAPACITY
182kg (28st)



WEIGHT
149kg (328lbs)



MAXIMUM
GRADIENT
10°



GROUND
CLEARANCE
100mm (3.9")

Leaflet Code: LL159

If you would like to view this product please contact your local stockist, details below.

