

# ROYALE 4 & 3

## SCOOTERS

*"The result is an exceptional driving experience"*



CE



DIGITAL LCD DISPLAY



ALLOY WHEELS



STYLISH DESIGN

The Drive Royale 3 and 4 are two of the finest mobility scooters available today, both models have an impressive performance; powered by 80amp batteries they have an excellent maximum range of 32 miles. These stylish scooters also feature a modern luminous LCD display, providing Speedo, High / Low turn speed, Power Indicator, Clock, Odometer, Temperature, Reverse Light and when necessary an Error Message will appear.

**drive**  
[www.drivemedical.co.uk](http://www.drivemedical.co.uk)

The Drive Royale Scooters have fully active independent front and rear suspension, extra-large 16" wheels with pneumatic tires and an ergonomic sport seat, the result is simple, an exceptional driving experience.

- Max range of 32 mile\*
- Top speed of 8mph\*
- 185kg (29st) weight limit
- All round suspension
- Digital dash
- Sliding and rotating seat
- Delta bar
- High level indicators.

- Stylish design
- 70ah batteries
- 4 colours available
- Automatic braking when steering.
- Large 16" diameter wheels
- Armrests are angle adjustable.

- Tiller adjustment by lever
- Twin rear view mirrors
- Road and pavement legal
- Comfortable seat with reclining facility
- Height adjustable headrest
- Optional rear pannier



## SPECIFICATION

Code: HW002  
Description: ROYALE 4




LENGTH  
158cm (62")



HEIGHT  
135cm (53")



WIDTH  
73cm (29")



MAX CAPACITY  
185kg (29st)



WEIGHT  
149kg (330lbs)



MAXIMUM GRADIENT  
?°




GROUND CLEARANCE  
150mm (6")

## SPECIFICATION

Code: HW001  
Description: ROYALE 3




LENGTH  
158cm (62")



HEIGHT  
135cm (53")



WIDTH  
73cm (29")




MAX CAPACITY  
185kg (29st)



WEIGHT  
149kg (330lbs)



MAXIMUM GRADIENT  
?°



GROUND CLEARANCE  
150mm (6")

### Colour Options



RED



WHITE



SILVER



BLACK

### Leaflet Code: LL158

If you would like to view this product please contact your local stockist, details across.

We have a continuous product improvement policy and consequently reserve the right to amend design and specification without prior notice. All sizes and weights stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate.