

PRISM SPORT

SCOOTER



"The high performance version of the popular Prism 4"

CE

"The Mercury range of Scooters continues to develop with the addition of the high performance Prism Sport. This Scooter is complete with exciting all new features, this includes; improved styling, a maximum range of 12 mile, a greater weight capacity and excellent new mud flaps."

drive
MEDICAL DESIGN & MANUFACTURING
www.drivemedical.co.uk

The Prism Sport is ideal for use indoors and outdoors where the terrain is forgiving (eg. smooth paths, shopping centres, etc). The Prism carries the Mercury brand, this ensures a focus on Comfort, Style, Performance and Reliability.



SPLITS EASILY FOR TRANSPORTATION & STORAGE

Superb Comfort, style & Performance...

- Maximum range of 12 miles on full battery charge*
 - Comfortable Bosun seat which swivels 360°
 - Excellent wireless connector means it can be split apart in seconds for easy transportation
 - Improved styling
 - Easy to remove battery pack fitted with larger 2 x 20ah sealed lead acid batteries
 - Excellent new mud flaps
 - New tiller head with headlight and LED battery gauge
 - Stylish and appealing body design
 - Larger 9" x 4" grey solid tyres providing greater ground clearance
 - Top speed of 4mph*
 - On board and in-house charger as standard
 - Intelligent electro magnetic brake system
 - Height adjustable seat
 - Removable seat, fold down tiller and integral lifting handles for easy transportation
 - Width and angle adjustable armrests
 - Front storage basket as standard
 - Anti tip wheels
 - Weight capacity 135kg (21st)
 - Delta bars available as optional extra
- *Dependent on user mass, terrain and battery condition.




INCREASED GROUND CLEARANCE

SPECIFICATION

Code: MS019
Description: PRISM SPORT




LENGTH
104cm (41")



HEIGHT
90cm (35.5")



WIDTH
55cm (21.5")



MAX CAPACITY
135kg (21St)



WEIGHT
47.2kg (106lbs)



MAXIMUM GRADIENT
8°



GROUND CLEARANCE
110mm (4")

Heaviest Component – 29.2kg (66lbs)



LARGER WHEELS



NEW 20AH BATTERY PACK



LED BATTERY GAUGE



NEW MUD FLAPS

Leaflet Code: LL143

If you would like to view this product please contact your local stockist, details below.