

# REGATTA

## SCOOTERS



*"The Regatta carries the Mercury brand: this ensures a focus on comfort, style, performance and reliability"*

The Regatta scooters have been developed for the more demanding journeys and include all the tried and tested benefits of the Mercury range. These easy to operate vehicles provide state of the art performance and will give years of dependable service. They have been developed with a focus on a stylish design, whilst also providing a quiet, smooth, comfortable and reliable ride every time.

CE

**drive**  
MEDICAL DESIGN & MANUFACTURING  
[www.drivemedical.co.uk](http://www.drivemedical.co.uk)



## REGATTA 6


4mph Pavement and 6mph maximum road speed settings\*  
–  
20 mile range\*  
–  
Weight carrying capacity 135kg (21 stone)

### SPECIFICATION

Code: MR6  
Description: REGATTA 6




LENGTH  
141cm



HEIGHT  
107cm



WIDTH  
64cm



MAX CAPACITY  
135kg (21 St)



WEIGHT  
96kg (211 lbs)



## REGATTA 8


4mph Pavement and 8mph maximum road speed settings\*  
–  
Weight carrying capacity 159kg (25 stone)  
–  
25 mile range\*

### SPECIFICATION

Code: MREG  
Description: REGATTA 8




LENGTH  
146cm



HEIGHT  
114cm



WIDTH  
67cm



MAX CAPACITY  
159kg (25 St)



WEIGHT  
117kg (257 lbs)




EASY TO USE PRECISION CONTROLS.



ONE TOUCH TILLER ADJUSTMENT



LARGE SLIDING & RECLINING SEAT



WIDTH & ANGLE ADJUSTABLE ARMRESTS



FULLY FUNCTIONING LIGHTS & INDICATORS FOR ROAD USE.



LARGE SEAT BACK STORAGE POCKET

## FEATURES

- Comfortable rotating Admirals seat with sliding and reclining facility.
- Easy to use controls.
- Adjustable seat position and seat back angle.
- Height adjustable headrest.
- Twin rear view mirrors.
- Anti slip carpet in footwell.
- Additional hand operated cable brake.
- Anti tip wheels.
- Width and angle adjustable armrests.
- One touch tiller adjustment.
- Road and pavement legal.
- Front and rear lights and indicators.
- Sturdy removable shopping basket.
- Large seat back storage pocket.
- Fully auto electric safe-braking system.
- Easy to use "Self Diagnostic" fault-finding system.
- Fully adjustable front and rear suspension.
- Delta Bar available as optional extra.
- Front and rear bumper.

\*Subject to conditions of use.

Leaflet Code: LL110 – If you would like to view this product please contact your local stockist, details right.



We have a continuous product improvement policy and consequently reserve the right to amend design and specification without prior notice. All sizes and weights stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate. \*Subject to Conditions of use.