

WEB: [WWW.SPORTRIDER.CO.UK](http://WWW.SPORTRIDER.CO.UK) EMAIL: [INFO@DRIVEMEDICAL.CO.UK](mailto:INFO@DRIVEMEDICAL.CO.UK)

# THE SPORT RIDER FROM DRIVE THE FUTURE HAS ARRIVED...

...A MILESTONE IN SCOOTER DESIGN, BUILT TO PUSH THE LIMITS.

**drive**

[WWW.DRIVEMEDICAL.CO.UK](http://WWW.DRIVEMEDICAL.CO.UK)



THE STYLISH NEW SPORT RIDER WAS INSPIRED BY SOME OF THE LEADING DESIGNS IN THE MOTORCYCLE INDUSTRY. THIS MOVES AWAY FROM THE TRADITIONAL TYPE OF SCOOTER, CREATING DESIRABILITY AND EXCITEMENT TO THE DRIVING EXPERIENCE.



## FEATURES

- EASY TO USE CONTROLS, TWIN DIAL DIGITAL DASH INCLUDING A TRIP COMPUTER WITH ODOMETER.
- FULLY AUTOMATIC EASY TO USE "SELF DIAGNOSTIC" FAULT-FINDING SYSTEM.
- ALL ROUND LIGHTING INCLUDING A LARGE HEADLIGHT, FRONT AND REAR INDICATORS AND REAR BRAKE LIGHTS.
- TELESCOPIC FRONT MOTORCYCLE SUSPENSION FOR OPTIMAL HANDLING.
- EXCELLENT MANOEUVRABILITY AND TURNING CIRCLE.
- COMFORTABLE ROTATING SPORT STYLE SEAT WITH SLIDING AND RECLINING FACILITY.
- WIDTH ADJUSTABLE ARMRESTS, HEADREST AND A LARGE SEAT BACK STORAGE POCKET.
- MOTORCYCLE STYLE WHEELS WITH CHROME MUD GUARDS AND ANTI TIP WHEELS.
- THROTTLE LEVER WITH AUTOMATIC SAFE-BRAKING SYSTEM PLUS AN ADDITIONAL HAND BRAKE.
- ANGLE ADJUSTABLE HANDLE BARS AND KEY IGNITION FOR SECURITY.
- REAR BUMPER PROVIDES ADDITIONAL PROTECTION.
- 4MPH PAVEMENT AND 8MPH MAXIMUM ROAD SPEED SETTINGS WITH A 31 MILE RANGE\*.
- LARGE ACCOMMODATING FOOTWELL WITH ANTI SLIP RUBBER MAT.

\*SUBJECT TO CONDITIONS OF USE




## SPECIFICATION

LENGTH: 165CM (65")  
WIDTH: 78CM (31")  
HEIGHT: 112CM (44")  
TURNING RADIUS: 157CM (62")  
GROUND CLEARANCE: 12CM (5")  
MAX RANGE: 31 MILES  
WEIGHT LIMIT: 180KG (28ST)  
MAX SPEED: 8MPH  
BATTERIES: 100AH 12V X 2  
BRAKE SYSTEM: ELECTROMAGNETIC

CLIMBING ANGLE: 12°  
CHARGING TIME: 12 HRS  
CONTROLLER: 160 AMP/DYNAMIC  
DISPLAY: DIGITAL  
COLOUR: BLACK & CHROME  
LEG ROOM: 51CM (20")  
SEAT TYPE: SPORTS STYLE  
SUSPENSION: FULL FRONT & REAR  
PANNIER: OPTIONAL EXTRA  
MIRRORS: LEFT & RIGHT AS STANDARD

FOR MORE INFORMATION PLEASE CONTACT YOUR LOCAL DEALER:



**drive**  
WWW.DRIVEMEDICAL.CO.UK

LEAFLET CODE: LL150