

# ST1

Scooter

Strider **ST#1**

*“Tight turning circle for ease of use indoors”*



CE

The ST1 Scooter is the lightest in the ST range and has an ultra compact design which is ideal for nipping to the shops, or for putting in the car boot to take on day trips and holidays. The streamlined ST1 offers everything you need for simple, independent convenience.

**drive**<sup>TM</sup>  
[www.drivemedical.co.uk](http://www.drivemedical.co.uk)

The STi's compact and stylish frame uniquely splits into two pieces, and can be further split into five smaller pieces for effortless lifting. Despite being the lightest scooter in the range, the STi has plenty of power to easily take you up slopes and hills, and is strong and stable for safety and security.



## Superb Comfort, Style & Performance...

- » Top speed of 4mph\*.
  - » Maximum range of 6 miles on a full battery charge\*.
  - » Height adjustable, hard wearing and fully cushioned profiled seat.
  - » Plug and wire free connections for quick and simple assembling and dismantling.
  - » Splits easily for storage and transportation.
  - » Width adjustable armrests and adjustable handlebars to suit comfort of user.
  - » Comfortable lightweight and rotating seat which swivels 360° for easy transfers on and off the scooter.
  - » Delta bars as standard.
  - » Removable seat, easy tiller adjustment and an integral lifting handle for easy transportation.
  - » Lightweight scooter weighing only 40kg (88lb).
  - » Horn to alert others of your presence for added safety.
  - » Fitted with reflectors and illuminated controls to improve outdoor visibility at dusk or when dark.
  - » Front removable storage basket as standard.
  - » Grey puncture proof non marking tyres.
  - » Rear anti tip wheels as standard.
  - » Anti roll back system safety device prevents scooter from rolling backwards on hills.
  - » Freewheel facility allows movement of scooter without switching the motor on.
  - » Easy to remove battery pack fitted with 2 x 12 Ah airline friendly batteries\*\*.
  - » Includes an on or off board charging facility as standard.
  - » Available in Ruby Red or Sapphire Blue.
  - » Optional extra's - Rear basket and additional storage options.
  - » Weight capacity 115kg (18st).
- \* Subject to conditions of use.  
\*\* Batteries included as standard.



Splits easily for storage and transportation

CODE: STiBE

DESCRIPTION: STi Scooter



Length:  
100cm (39")



Height:  
89cm (35")



Width:  
51cm (20")



Max Capacity:  
115kg (18st)



Weight:  
40kg (88lb)



Max Gradient:  
6°

### Colour Options



Sapphire Blue

Ruby Red

Leaflet Code: LL191

If you would like to view this product please contact your local stockist, details below.



Delta bars as standard