

# ST6

Scooter

Strider **ST6**

*“Excellent on and off road!”*



CE

The ST6 Scooter offers a modern design and a high performance allowing the scooter to be used both on and off road. The scooter is fitted with a range of safety features, including two rear view mirrors and LED lights.

**drive**<sup>TM</sup>  
[www.drivemedical.co.uk](http://www.drivemedical.co.uk)

The ST6 is a perfectly formed scooter with its sweeping curves and enviable design, and offers an excellent weight capacity of 31st and a maximised range of 31 miles. The scooter features extra ground clearance and oil damped suspension, both of which make your journey smoother and more comfortable. On top of all these features, this scooter has an excellent grade climb to ensure you can effortlessly go up hills.



- » Top speed of 8mph\*.
- » Maximum range of 31 miles on a full battery charge\*.
- » Offers superb comfort with a luxurious, reclining, angle adjustable and sliding seat with a padded adjustable headrest and lumbar support.
- » Luxurious Captains seat with two tone cover.
- » High / low speed selector to switch between 4mph (pavement only) and 8mph (road legal).
- » High ground clearance of 10cm (5") for a smooth ride over rough terrain.
- » Fully functioning indicators and adjustable twin rear view mirrors for extra visibility.
- » Oil damped suspension ensures a smooth and comfortable ride.
- » Backlit LCD illuminated controls.
- » Fitted with innovative Rhino 2 control system for a smooth, powerful ride.
- » Comfortable lightweight and rotating seat which swivels 360° for easy transfers on and off the scooter.
- » Removable seat and pneumatic tiller adjustment for easy transportation.
- » Comfortable height adjustable seat with adjustable armrests.
- » Delta bars as standard.
- » Four tie down points for safe transportation.
- » Horn to alert others of your presence for added safety.
- » Fitted with reflectors and LED front, rear and brake lights to improve outdoor visibility at dusk or when dark.
- » Power / Self-Diagnostic and Hazard warning lights for added safety.
- » Alloy wheels with black pneumatic tyres for a smoother ride.
- » Rear anti tip wheels as standard
- » Anti roll back system safety device prevents scooter from rolling backwards on hills.
- » Freewheel facility allows movement of scooter without switching the motor on.
- » Easy to remove 2 x 75 Ah airline friendly batteries\*\*.
- » Includes an on board charging facility as standard.
- » Emergency brake for added safety.
- » Available in Silver.
- » Additional storage in the modesty panel and seat pocket.
- » Wide range of optional extra's including rear basket and additional storage options.
- » Weight capacity 197kg (31st).

\* Subject to conditions of use.  
\*\* Batteries included as standard.

CODE: **ST6SBE**  
DESCRIPTION: **ST6 Scooter**



Backlit LCD illuminated controls



Length:  
160cm (63")




Height:  
128cm (50")



Width:  
72cm (28")



Max Capacity:  
197kg (31st)



Ground Clearance:  
10cm (4")



Max Gradient:  
9°



Alloy wheels with black pneumatic tyres for a smoother ride

#### Colour Options



Silver

#### Leaflet Code: LL197

If you would like to view this product please contact your local stockist, details below.