

# Alcora

Model No. **S740**

**ROMA**



## Special Features

- Detachable battery pack for off board charging.
- Quickly dismantles without tools for easy transportation and storage.
- Non marking puncture proof tyres.
- Adjustable height seat.
- Flip up and width adjustable armrests.
- Adjustable angle tiller.



## Specifications



Model No.*	<b>S740</b>	Maximum User Weight	<b>100 kgs / 16 stone</b>
Class (Invalid Carriage)	<b>2</b>	Battery Size 12v x 2	<b>12 Ah x 2 (10 A/h)</b>
Overall Length	<b>955 mm / 37.5"</b>	Range**	<b>10 km / 6 miles</b>
Overall Width	<b>545 mm / 21.5"</b>	Weight	<b>35 Kgs / 77 lbs</b>
Overall Height	<b>920 mm / 36"</b>	Weight of Heaviest Part (Chassis Assembly)	<b>21.5 kgs / 47 lbs</b>
Seat Width	<b>410 mm / 16"</b>	Weight of Battery Pack	<b>7 kgs / 15 lbs</b>
Seat Depth	<b>340 mm / 13.5"</b>	Weight of Seat	<b>6.5 kgs / 14 lbs</b>
Seat Height (to ground)	<b>550mm - 600mm / 22"-24"</b>	Maximum Speed	<b>4 mph / 6kmph</b>
Armrest Height	<b>190 mm / 7.5"</b>	Wheel Size (Front / Back)	<b>180 mm / 7" 200 mm / 8"</b>
Backrest Height (from Seat)	<b>320 mm / 12.5"</b>	Turning Radius	<b>1105 mm / 43.5"</b>
Seat Height (to Floorpan)	<b>430 - 480mm / 17"-19"</b>	Battery Charger	<b>1.8 Amp</b>
Safe Working Slope	<b>6° / 10%</b>	Ground Clearance***	<b>90mm / 3.5"</b>

\* Lower limb amputees should refer to user manual before use. \*\* Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel. \*\*\* Measurement to underside of chassis